

# NOR-LAKE SCIENTIFIC

## GENERAL PURPOSE LABORATORY REFRIGERATOR

### UPRIGHT AUTO DEFROST

#### Standard Product Features:

- Operating temperature: +4°C Refrigerator
- Auto Defrost, Forced Air
- Adjustable Temperature Control
- 2 Adjustable Wire Shelves
- 2 Fixed Wire Shelves
- 5 Full-Width Fixed Door Bins (1 Tall Storage)
- Bottom Shelf
- Interior Light
- 2 Adjustable Leg Levelers
- Automatic Door Closer
- White Interior, Exterior, Door and Handle
- Complete Cabinet is Foamed-In-Place with CFC Free High Density Polyurethane Foam Insulation
- Warranties: 12 Months Parts



**Digital Thermometer Alarm (optional)**

#### Optional Features:

- Digital Thermometer Monitor with Alarm (Field Install)



**LR161WWW/0**



**LR161WWW/0**



NOR-LAKE SCIENTIFIC  
A Division of  
NOR-LAKE, INCORPORATED  
727 Second Street  
P.O. Box 248  
Hudson, Wisconsin 54016  
800-477-5253  
715-386-2323  
800-388-5253 PARTS  
715-386-4290 FAX  
[www.norlake.com](http://www.norlake.com)

